



A BURN
TAKES A SECOND
TO OCCUR
& A LIFETIME TO
HEAL IF TREATED
INCORRECTLY



CATALOGUE
2022



MDSAP



TGA
Health Safety
Regulation

SAHPRA
Health Practice
Regulatory Authority

WORLD CLASS FACTORY IN JOHANNESBURG, SOUTH AFRICA

We are a South African based medical manufacturer of Burnshield Emergency Burn Care and First Aid Kits, audited annually with CE & ISO accreditation.

Levtrade International (PTY) LTD was founded in 1990 with a single product line, Burnshield®. Burnshield® is an international go-to emergency burn care range of premium products that provide an immediate and effective method of treating burn injuries in the acute (emergency) phase.

Burnshield® is manufactured and distributed from our world class factory in Johannesburg, South Africa which meets the requirements of and holds certificates of approval to market. These accreditations have afforded our company the opportunity to export into more than 42 countries worldwide.

Burnshield® is endorsed and accepted internationally through independent laboratory tests, research and acceptance in key Hospital & EMS protocols. This paved the way for the development into various markets including: Medical, Consumer, Industrial, Military, Police/Traffic, Correctional Services, First Aid Kit Manufacturers, Hospital and EMS/Fire departments. Training & Educational items have been introduced to compliment the range.



As the company grew and opened up new markets, the demand for a quality First Aid & Surgical Range became a priority and a perfect compliment to the Burnshield range, Levtrade designed and produced an extensive range of First Aid & Trauma Kits to cover every purpose including Domestic, Factories, Hotels, Sport, Retail, Medical Professionals and Corporates.

The product range expanded to include Solace® Suncare products endorsed by CANSA South Africa for UVA and UVB protection. Solace® Suncare is the proud sponsor of IRONMAN South Africa.

Our belief in providing quality at all times afforded us the opportunity to distribute our own brand of Sport Strapping, Tapes and Kits under the PROSTRAP® brand. This has been well received by Sport Doctors, Sports Teams, Physio's, Clubs and Schools.

Our manufacturing and warehouse facilities ensure we provide quality and service for a variety of Professional, Occupational Health, Corporate, Retail and Relief organisations.

We offer personal service, premium support and delivery to our clients worldwide. It is due to the continued support of our many loyal clients and the team effort required in-house that allows us to enjoy the success we do both locally and internationally.

Thank you for taking the time to browse our catalogue and we look forward to working with you in the future.

BURNSHIELD® products meet the requirements of and hold certificates to: ISO 13485 and CE (European mark for the highest standard of manufacture and quality), MDSAP (Canada), FDA (USA), TGA (Australia), and SAHPRA (South African Health Products Regulatory Authority).




MDSAP
MEDICAL DEVICE SINGLE AUDIT PROGRAM



TGA
Health Safety
Regulation


SAHPRA
South African
Health Products
Regulatory Authority



KEY BENEFITS

- STERILE
- MINIMISES TRAUMA & SKIN DAMAGE
- MOISTENS, COOLS & SOOTHES THE WOUND
- PHYSICAL PROTECTION AGAINST INFECTION
- NON-ADHERENT
- NON TOXIC
- NON-IRRITANT
- SAFE TO USE ON CHILDREN
- PRODUCTS ARE STERILIZED BY GAMMA IRRADIATION

BURNSHIELD® IS SUITABLE FOR ALL TYPES OF BURNS



HYDROGELS



HYDROGELS

SUITABLE FOR

Sunburn
Minor burns and scalds
Difficult to reach areas: e.g. behind the ears,
creases of the face and groin area
Rehydrate the Burnshield Dressing



125 ml
SPRAY BOTTLE
4,50 oz



75 ml
SPRAY BOTTLE
2,70 oz



50 ml
SQUEEZE BOTTLE
1,80 oz



50 ml
SPRAY BOTTLE
1,80 oz



25 ml
TUBE
0,85 oz



10 ml
TUBE
0,35 oz



3,5 ml SACHET
1/8 oz



0,9 g SACHET
1/32 oz



CARTRIDGES
1.2 kg 42,32 oz
150 g 5,29 oz

SUITABLE FOR
EMS
Military
Industry

HYDROGEL APPLICATION

- * Ensure to extend application beyond the burn area
- ** Re-apply to prevent the wound from dehydrating
- *** Seek medical attention



Apply the Burnshield Hydrogel thoroughly



Secure dressing with non-absorbant bandage



KEY BENEFITS

Compact. Light weight.
Sterile.
Available in a variety of different packaging: bottles, tubes and sachets.
Items are individually labelled with instructions for use.

Both the Burnshield Dressing and Burnshield Hydrogel can be used independently or in conjunction with each other

ITEM	SINGLE UNIT					SHIPPER				CARTON				
	REF	BARCODE	MIN QTY	W x D x H (mm)	(g)	BARCODE	UNITS	W x D x H (mm)	(g)	BARCODE	UNITS	W x D x H (mm)	(kg)	QTY PER PALLET SEA FREIGHT (1 x 1 x 1.8m (max height))
125 ml / 4,5 oz Spray Bottle	550003	6007645000192	48	45 x 35 x 170	153	16007645000199	6			36007645000193	48	300 x 300 x 190	8,18	72
75 ml / 2,7 oz Spray Bottle	550082	6007645004305	48	53 x 35 x 112	94	26007645004309	6			16007645004302	48	340 x 310 x 130	5,14	108
50 ml / 1,8 oz SPRAY BOTTLE	550020	6007645006040	48	40 x 25 x 130	71			85 x 90 x 110	399	16007645006047	48	250 x 240 x 160	3,7	160
50 ml / 1,8 oz Squeeze Bottle	550000	6007645000123	72	40 x 25 x 108	65	16007645000120	6			26007645000127	72	240 x 170 x 230	5	168
25 ml / 0,85 oz TUBE	550009	6007645005302	72	45 x 25 x 105	31	36007645005303	6	63 x 83 x 106	193	46007645005300	72	360 x 240 x 130	2,62	120
10 ml / 0,35 oz TUBE	550010	6007645004596	500	35 x 22 x 70	775	26007645004590	50		20	36007645004597	500	390 x 330 x 280	8,07	40
3,5 ml / 1/8 oz SACHET	550004	6007645000215	1200	45 x 2 x 85	5	-	-	-	-	16007645000212	1200	415 x 300 x 210	6,2	64
0,9 g / 1/32 oz SACHET	550021	6007645006057	2400	45 x 2 x 60	2					16007645006054	2400	320 x 285 x 210		64
1.2 kg 42,32 oz	550022	6007645007078	6			16007645007075	6							
150 g 5,29 oz	550023	6007645007153	30	50 x 35 x 115	90	16007645007150	30							

THE BURNSHIELD DRESSING IS THE PRIMARY PRODUCT OF CHOICE:

A fibreless elastic open-cell foam sheet dressing soaked in Burnshield sterile hydrogel, supplied in a foil pouch.


KEY BENEFITS

- Optimal heat dissipation
- Natural wound covering
- Keeps the burn area protected and moist
- Conforms to wound area
- Gel retention
- Limits bacteria exposure to the wound

Both the Burnshield Dressing and Burnshield Hydrogel can be used independently or in conjunction with each other



BURNSHIELD DRESSING APPLICATION



Open the sterile package

Remove dressing from package

Apply dressing over the wound

Secure dressing with non-absorbant bandage

* Re-apply Burnshield Hydrogel to prevent the wound from dehydrating.
** Seek medical attention.



100 mm x 100 mm
4" x 4"

200 mm x 200 mm
8" x 8"

200 mm x 450 mm
8" x 18"

600 mm x 400 mm
24" x 16"

200 mm X 450 mm
FACE MASK (PRE-CUT)
8" x 18"

600 mm x 400 mm
FACE MASK (PRE-CUT)
24" x 16"

50 mm x 1 m LIMB
2" x 40"

25 mm x 500 mm DIGIT
1" x 20"

ITEM	SUITABLE FOR	SINGLE UNIT					SHIPPER					CARTON				
		REF	BARCODE	MIN QTY	W x D x H (mm)	(g)	BARCODE	UNITS	W x D x H (mm)	(g)	BARCODE	UNITS	W x D x H (mm)	(kg)	QTY PER PALLET SEA FREIGHT 1x1.8m (max height)	
100 mm x 100 mm 4" x 4"	• Small Sized Burns • Fingers	901111	6007645001380	250	118 x 3 x 123	47	16007645001387	25	120 x 210 x 153	1279	36007645001381	250	350 x 425 x 325	13,305	30	
200 mm x 200 mm 8" x 8"	• Medium Sized Burns • Facial Burns	900906	6007645001397	50	230 x 3 x 198	142	16007645001394	10	35 x 245 x 205	1542	26007645001391	50	230 x 255 x 285	7,96	80	
200 mm x 450 mm 8" x 18"	• Medium Sized Burns • Facial Burns	901001	6007645001410	30	250 x 5 x 230	293	16007645001417	5	35 x 250 x 225	1698	26007645001414	30	230 x 255 x 285	9,53	80	
600 mm x 400 mm 24" x 16"	• Large Sized Burns • Facial Burns	901007	6007645001403	12	435 x 5 x 300	568	-	-	-	-	16007645001400	12	325 x 480 x 950	7,25	108	
50 mm x 1 m LIMB 2" x 40" SINGLE	• Limb & neck • Circumferential Burns	901131	6007645001427	50	250 x 15 x 130	162	16007645001424	10	65 x 260 x 210	1645	26007645001421	50	230 x 255 x 285	8,84	80	
25 mm x 500 mm DIGIT 1" x 20"	• Fingers & Toes • Circumferential Burns	901130	6007645001434	200	120 x 5 x 94	45	16007645001431	20	125 x 122 x 127	966	26007645001438	200	325 x 385 x 290	10,495	40	
FACE MASK 200 mm x 450 mm (PRE-CUT) 2" x 40"	• Facial Burns	901000	6007645001458	30	250 x 5 x 230	301	16007645001455	5	35 x 250 x 225	1819	26007645001452	30	230 x 255 x 285	9,44	80	
FACE MASK 600 mm x 400 mm (PRE-CUT) 2" x 40"	• Facial Burns	901135	6007645003957	12	435 x 5 x 300	560	-	-	-	-	16007645003954	12	325 x 480 x 950	7,355	108	



FIRE TRAUMA BLANKET

WOOL



Burnshield's original ground-breaking product, Fire Trauma Blanket was developed in conjunction with the Wool Board and the CSIR. The Fire Trauma Blanket is used as a large area covering, escape from a burning area and to extinguish flames on a person or object. It provides pain relief and healing while victims are being transported to emergency care units.



50% of all burn injuries occur in the home, and the vast majority of the victims are young children. This unique product is effective in treating severe upper body burns on children from 6 months to 5 years old. This Cotton Burn Vest is manufactured from fibreless cotton for comfort and maximum gel retention. The Burnshield Burn Vest has been tested and is supported by the Paediatric Burns Unit at the Chris Hani Baragwanath Hospital.



FIBRELESS COTTON

PAEDIATRIC BURN VESTS



ITEM	CODE	BARCODE	CONTENTS	W x D x H (mm)	MIN QTY	QTY PER PALLET SEA FREIGHT 1x1x1.8m (max height)	KEY FEATURES
1.2 m x 1.6 m 48" x 64" SINGLE	881003	6007645000222	x1 Blanket 1.2 m x 1.6 m	500 x 185 x 185	1	50	<ul style="list-style-type: none"> Wool base Wound coverage Heat protection Supplied in PVC barrel shoulder bag Suitable for: Industry & EMS Military Mining
1.2 m x 1.6 m 48" x 64" TWIN	881002	6007645002073	x2 Blanket 1.2 m x 1.6 m	500 x 185 x 185	1	50	
2.45 m x 1.6 m 96" x 64"	881001	6007645000161	x1 Blanket 2.45 m x 1.6 m	500 x 185 x 185	1	50	
1.6 m x 1.8 m 64" x 71"	881010	6007645007917	x1 Blanket 1.6 m x 1.8 m	500 x 185 x 185	1	50	

ITEM	SUITABLE FOR	SINGLE UNIT				CARTON			
		REF	BARCODE	W x D x H (mm)	(g)	BARCODE	QTY	W x D x H (mm)	(Kg)
6 - 12 MONTHS PRE-CUT	<ul style="list-style-type: none"> Torso, Back & Arms Large Area Burns 6 - 12 Months 	881080	6007645004329	435 x 15 x 320		16007645004326	5		
24 - 36 MONTHS PRE-CUT	<ul style="list-style-type: none"> Torso, Back & Arms Large Area Burns 2 years - 3 Years 	881085	6007645004336	435 x 15 x 320		16007645004333	5		
4 - 5 YEARS PRE-CUT	<ul style="list-style-type: none"> Torso, Back & Arms Large Area Burns 4 years - 5 Years 	881090	6007645004343	435 x 15 x 320		16007645004340	5		

The flexible fibreless cotton based Contour Burn Dressing was developed as an extension to the professional range as a half body and full body dressing, which makes this cost effective and unique for the hospital, EMS, military and occupational health.



ITEM	CONTENTS	SINGLE UNIT					SHIPPER			CARTON				
		REF	BARCODE	MIN QTY	W x D x H (mm)	(g)	BARCODE	W x D x H (mm)	(g)	BARCODE	UNITS	W x D x H (mm)	(Kg)	QTY PER PALLET SEA FREIGHT 1x1x1.8m (max height)
1 m x 1 m 40" x 40" HALF BODY	x1 Dressing 1 m x 1 m	881006	6007645003520	5	430 x 40 x 350	1231	6007645003520			16007645003527	5	320 x 260 x 150	6,4	132
1 m x 2 m 40" x 80" FULL BODY	x1 Dressing 1 m x 2 m	881007	6007645003537	5	430 x 40 x 350	2256				16007645003534	5	320 x 270 x 230	12,1	63
1 m x 1 m 40" x 40" IN A CYLINDRICAL BAG	x1 Dressing 1 m x 1 m x1 Cylindrical Bag	881060	6007645007641	5	340 x 110 x 110			340 x 115 x 120	1333	16007645007648	5	370 x 255 x 290	7,5	50
1 m x 2 m 40" x 80" IN A NYLON BAG	x1 Dressing 1 m x 2 m x1 Nylon Bag	881070	6007645007658	5	BAG: 265 x 90 x 450			285 x 55 x 350	2684	16007645007655	5	370 x 255 x 290	13,8	50
BURN KIT 1 m x 2 m 40" x 80" IN A NYLON BAG	6 ITEMS x2 Hydrogel 125 ml x1 Dressing 1 m x 2 m x1 Cohesive Bandage 50 mm x1 Cohesive Bandage 75 mm x1 Nylon Bag	881008	6007645007665	5	BAG: 265 x 90 x 450					16007645007662	5			
BURN KIT 1 m x 1 m 40" x 40" DOUBLE IN A NYLON BAG	7 ITEMS x2 Hydrogel 125 ml x2 Dressing 1 m x 1 m x1 Cohesive Bandage 50 mm x1 Cohesive Bandage 75 mm x1 Nylon Bag	881009	6007645007672	5	265 x 90 x 450					16007645007679	5			

FIBRELESS COTTON

CONTOUR DRESSINGS





RETAIL

RETAIL

Burnshield products are compact and light weight. They are packaged for intended hanging and shelf display in-store. All header card items contain a Euro-slot Hanging Tag which easily hang the products for display in:

- Pharmacy
- Retail
- FMCG



MINI BURN KIT SINGLE



MINI BURN KIT DISPENSER OF 6'S



HOME BURN KIT 10B



75 ml SPRAY BOTTLE + SOFT 1 BLISTER PACK



50 ml SQUEEZE BOTTLE + SOFT 1 BLISTER PACK



OWN BRANDING OFFER

With the minimum order quantities being easily achievable for many of our products, this service can be a perfect option to promote your own brand. We manufacture and design to your requirements with our in-house design department. We guarantee fast and efficient service.



25 ml TUBE BLISTER PACK OF 2'S



25 ml TUBE DISPENSER OF 12'S



3,5 ml SACHET BLISTER PACK OF 5'S



3,5 ml SACHET BOX OF 10'S



100 mm x 100 mm TWINPACK



100 mm x 100 mm DISPENSER OF 6'S



100 mm x 100 mm DISPENSER OF 15'S

ITEM	CONTENTS	SINGLE UNIT					SHIPPER				CARTON				
		REF	BARCODE	MIN QTY	W x D x H (mm)	(g)	BARCODE	UNITS	W x D x H (mm)	(g)	BARCODE	UNITS	W x D x H (mm)	(kg)	QTY PER PALLET SEA FREIGHT 1x1x1.8m (max height)
MINI BURN KIT SINGLE	x1 Dressing 100 mm x 100 mm x1 Cohesive Soft 1 Foam Bandage x2 Hydrogel Sachet 3,5 ml	550013	6007645001540	100	120 x 25 x 150	92			120 x 25 x 153	74	36007645001541	100	520 x 340 x 330	1,09	
MINI BURN KIT DISPENSER OF 6'S	x6 Mini Burn Kits	550014	6007645002035	16	OPEN: 122 x 157 x 230 CLOSED: 122 x 157 x 155	572			125 x 155 x 130		26007645001544	16			
100 mm x 100 mm TWINPACK	x2 Dressing 100 mm x 100 mm	550060	6007645004527	100	120 x 17 x 155	110					16007645004524	100			
100 mm x 100 mm DISPENSER OF 6'S	x6 Dressing 100 mm x 100 mm														
100 mm x 100 mm DISPENSER OF 15'S	x15 Dressing 100 mm x 100 mm	550006	6007645000376	12	OPEN: 120 x 118 x 210 CLOSED: 120 x 118 x 150	740	-		125 x 110 x 125	738	56007645001385	12	385 x 230 x 285	7,9	32
HOME BURN KIT 10B	x1 Dressing 100 mm x 100 mm x4 Hydrogel Sachet 3,5 ml x1 Hydrogel 50 ml	550008	6007645001564	15	112 x 45 x 225	168					26007645001568	15			
75 ml SPRAY BOTTLE + SOFT 1 BLISTER PACK	x1 Hydrogel 75 ml x2 Hydrogel Sachet 3,5 ml x1 Soft Cohesive Bandage	550075	6007645004558	50	110 x 50 x 210										
50 ml SQUEEZE BOTTLE + SOFT 1 BLISTER PACK	x1 Hydrogel 50 ml x1 Soft Cohesive Bandage	550080	6007645004565	50	110 x 40 x 207										
25 ml TUBE BLISTER PACK OF 2'S	x2 Hydrogel 25 ml	550018	6007645005753	50	210 x 110 x 30	90	16007645005750	5			26007645005757	50			
25 ml TUBE DISPENSER OF 12'S	x12 Hydrogel 25 ml	550016	16007645005309	12	OPEN: 135 x 135 x 173 CLOSED: 135 x 135 x 106				135 x 135 x 116	417	26007645005306	12	380 X 330 X 290	5,74	
3,5 ml SACHET BLISTER PACK OF 5'S	x5 Hydrogel Sachet 3,5 ml	550065	6007645004534	100	85 x 20 x 170	36	16007645004531	10	85 x 25 x 170	36	36007645004535	100	410 x 290 x 210	6,03	
3,5 ml SACHET BOX OF 10'S	x10 Hydrogel Sachet 3,5 ml	550005	6007645000017	50	48 x 46 x 85	50					36007645000018	50	480 x 325 x 95	3	108
3,5 ml SACHET BOX OF 25'S	x25 Hydrogel Sachet 3,5 ml	550017	6007645005746	40	48 x 95 x 85	125	-	-	-	-	26007645005740	40	260 x 230 x 370	5,56	50
0,9 g SACHET DISPENSER OF 144'S	x144 Hydrogel Sachet 0.9 g	550090	26007645006051	16	139 x 70 x 171	297					36007645006058	16			



3,5 ml SACHET BOX OF 25'S



0,9 g SACHET DISPENSER OF 144'S

SUITABLE FOR:

- Home Care
- Travel
- Sun Burn
- Rehydration
- Small + Medium Burns
- First Aid Kits
- Consumer
- Retail
- Fast Food Outlets

MUNDIAL

PREMIUM RANGE

MUNDIAL



PREMIUM RANGE

Burnshield Mundial is a premium packaged range of Burnshield products intended for retail markets. soft foil packs are compact and light weight. Items are individually labelled with instructions for use.



125 ml SPRAY BOTTLE
4,50 oz



50 ml SQUEEZE BOTTLE
1,80 oz



25 ml TUBE
0,85 oz



TRAVEL PACK



HOME / TRAVEL BURN KIT



TWO - PACK



BURN KIT



EMERGENCY BURN PACK

ITEM	REF	BARCODE	CONTENTS	MIN QTY	SINGLE UNIT		SHIPPER				CARTON			
					W x D x H (mm)	(g)	BARCODE	UNITS	W x D x H (mm)	(kg)	QTY PER PALLET	SEA FREIGHT 1x1x1.8m (max height)		
125 ml 4,50 oz SPRAY BOTTLE	550024	6007645007368	x1 Hydrogel 125 ml Spray Bottle	48	45 x 35 x 170	153				16007645007365	48	300 x 300 x 190	8,18	72
75 ml 2,7 oz SPRAY BOTTLE	550031	6007645008266	x1 Hydrogel 75 ml Spray Bottle	48	54 x 32 x 110	90				16007645008263	48	340 x 310 x 130	5,14	108
50 ml 1,80 oz SQUEEZE BOTTLE	550026	6007645007375	x1 Hydrogel 50 ml Squeeze Bottle	72	40 x 25 x 115	65	16007645007372	6		26007645007379	72	240 x 170 x 230	5,035	168
25 ml 0,85 oz TUBE	550025	6007645007054	x1 Hydrogel 25 ml Tube	72	45 x 25 x 105	31	16007645007051	6		26007645007058	72	360 x 240 x 130	2,62	120
TRAVEL PACK	550005F	6007645006361	x10 Hydrogel Sachet 3,5 ml	80	108 x 30 x 150	50	16007645006368	5		26007645006365	80			
HOME / TRAVEL BURN KIT	550008F	6007645006354	x1 Dressing 100 mm x 100 mm x4 Hydrogel Sachet 3,5 ml x1 Hydrogel 50 ml x1 Conform Bandage 50 mm x 2 m	80	138 x 25 x 215		16007645006351	5		26007645006358	80			
TWO-PACK	550060F	6007645006347	x2 Dressing 100 mm x 100 mm	80	138 x 25 x 215	110	16007645006344	5		26007645006341	16			
BURN KIT	550013F	6007645006330	x1 Dressing 100 mm x 100 mm x2 Hydrogel Sachet 3,5ml x1 Conform Bandage 50 mm x 2 m	80	138 x 25 x 215	70	16007645006337	5		16007645006334	80			
EMERGENCY BURN PACK	900802F	6007645006378	x2 Dressing 100 mm x 100 mm x2 Dressing 200 mm x 200 mm x2 Dressing 600 mm x 400 mm	5	290 x 35 x 400					16007645006375	5			





32 items

+ Back Straps

SUITABLE FOR
Occupational Health, EMS / Fire
Mountain Rescue
Military
Mining
Rescue Services

CONFINED RESCUE KIT

REF	900815
BARCODE	6007645005609
SIZE IN MM	270 x 140 x 500
DETAIL	Nylon back pack with back straps designed for search and rescue access into confined areas. Velcro attachment to secure a Trauma Blanket / Contour Dressing. Back flap for easy access. 3 x pull-out inner pouch system. Items individually labelled with instructions for use.
CONTENTS: 32 ITEMS	x8 Burnshield Dressing 100 mm x 100 mm x3 Burnshield Dressing 200 mm x 200 mm x2 Burnshield Dressing 600 mm x 400 mm x3 Burnshield Digit Dressing 25 mm x 0,5 m x2 Burnshield Limb Dressing 50 mm x 1 m x3 Burnshield Hydrogel Spray 125 ml x1 Burnshield Contour 1 m x 1 m x2 Crepe Bandage 150 mm x1 Plastipore Plastic Tape 25 mm x3 Cohesive Soft 1 Foam Bandage x2 Pair of Gloves x1 First Aid Scissors (Bandage) x1 Rescue Blanket Adult
QTY / CASE	1



45 items

+ D Handles

SUITABLE FOR
Industry & EMS
Petro-Chemical
Military
Mining
Disaster

COMBI KIT

REF	900800
BARCODE	6007645002325
SIZE IN MM	540 x 120 x 395
DETAIL	Sturdy rigid nylon case with D-handles, suitable for major burn incidents. Items individually labelled with instructions for use.
CONTENTS: 45 ITEM	x5 Burnshield Dressing 100 mm x 100 mm x5 Burnshield Dressing 200 mm x 200 mm x3 Burnshield Dressing 600 mm x 400 mm x5 Burnshield Digit Dressing 25 mm x 0,5 m x3 Burnshield Limb Dressing 50 mm x 1 m x4 Burnshield Hydrogel Spray 125 ml x2 Burnshield Cartridge 1.2 kg x1 Burnshield Contour 1 m x 1 m x1 Burnshield Contour 1 m x 2 m x2 Burnshield Cartridge 150 g x2 Crepe Bandage 100 mm x1 Plastipore Plastic Tape x4 Cohesive Soft 1 Foam Bandage x2 Pair of Gloves x1 First Aid Scissors (Bandage) x1 Rescue Blanket Adult x1 Syringe 20 ml x2 Purified Water - 1 Litre
QTY / CASE	1



16 items

+ D Handles

SUITABLE FOR
Home / Travel
Retail
Industry

RESCUE KIT

REF	900813
BARCODE	6007645002028
SIZE IN MM	145 x 115 x 320
DETAIL	Nylon bag. Portable with D-handles, professional and comprehensive. Items individually labelled with instructions for use.
CONTENTS: 16 ITEMS	x2 Burnshield Dressing 100 mm x 100 mm x2 Burnshield Dressing 200 mm x 200 mm x2 Burnshield Dressing 600mm x 400 mm x1 Burnshield Digit Dressing 25 mm x 0,5 m x1 Burnshield Limb Dressing 50 mm x 1 m x1 Burnshield Hydrogel Spray 125 ml x2 Conform Bandage 50 mm x1 Plastipore Plastic Tape x1 Soft 1 Cohesive Bandage x1 Pair of Gloves x1 First Aid Scissors (Bandage) x1 Rescue Blanket Adult
QTY / CASE	5



20 items

+ D Handles

SUITABLE FOR
EMS
Military
Industry

EMERGENCY KIT

REF	900817
BARCODE	6007645007924
SIZE IN MM	140 x 70 x 240
DETAIL	Nylon bag. Portable with D-handles, professional and comprehensive. Items individually labelled with instructions for use.
CONTENTS: 20 ITEMS	x2 Burnshield Dressing 100 mm x 100 mm x1 Burnshield Dressing 200mm x 200 mm x3 Burnshield Digit Dressing 25 mm x 0,5 m x1 Burnshield Limb Dressing 50 mm x 1 m x10 Burnshield Hydrogel Sachet 3,5 ml x1 Burnshield Hydrogel Spray 125 ml x1 Plastipore Plastic Tape x1 Soft 1 Cohesive Bandage
QTY / CASE	5



33 items

+ D Handles

SUITABLE FOR
EMS
Military
Industry

DISASTER KIT

REF. CODE	900700
BARCODE	6007645008051
SIZE IN MM	
DETAIL	Portable with D-handles. Professional and comprehensive.
CONTENTS: 11 ITEMS	x8 Burnshield Cartridge 1.2 kg x3 Burnshield Cartridge 150 g
QTY / CASE	1
REPLENISHMENT STOCK	
BARCODE	16007645002056
QTY / CASE	4



33 items

+ D Handles

SUITABLE FOR
EMS
Military
Industry

RESPONDER KIT

CODE	900814
SIZE IN MM	500 x 120 x 250
DETAIL	Portable with D-handles. Professional and comprehensive. Nylon Bag with pull-out pouches. Items individually labelled with instructions for use.
CONTENTS: 33 ITEMS	x8 Burnshield Dressing 100 mm x 100 mm x3 Burnshield Dressing 200 mm x 200 mm x2 Burnshield Dressing 600 mm x 400 mm x3 Burnshield Digit Dressing 25 mm x 0,5 m x2 Burnshield Limb Dressing 50 mm x 1 m x3 Burnshield Hydrogel Spray 125ml x2 Conform Bandage 75 mm x1 Plastipore Plastic Tape x4 Soft 1 Cohesive Bandage x2 Pair of Gloves x1 First Aid Scissors (Metal) x1 Rescue Blanket Adult x1 Syringe 20 ml

QTY / CASE 4



12 items

+ D Handles

SUITABLE FOR
Restaurants / Canteen
Hospitality
Industry
Lodges

CATER KIT

CODE	903001
SIZE IN MM	255 x 85 x 300
DETAIL	Catering / hospitality industry burn kit. Nylon bag. Portable with D-handles. Items individually labelled with instructions for use.
CONTENTS: 12 ITEMS	x2 Burnshield Dressing 100 mm x 100 mm x2 Burnshield Dressing 200 mm x 200 mm x1 Burnshield Hydrogel 125 ml x1 Burnshield Contour Dressing 1 m x 1 m x2 Conform Bandage 50 mm x1 Plastipore Plastic Tape x1 Pair of Gloves x1 First Aid Scissors (Bandage) x1 Rescue Blanket Adult

QTY / CASE 5



4 items

PVC Folder

SUITABLE FOR
EMS / Fire
Trauma Crash-Cart
Industry

MULTI PACK

CODE	900802
SIZE IN MM	375 x 40 x 300
DETAIL	Variety of dressings to suit most burn wounds. PVC folder for easy access and storage. Items individually labelled with instructions for use.
CONTENTS: 6 ITEMS	x2 Burnshield Dressing 100 mm x 100 mm x2 Burnshield Dressing 200 mm x 200 mm x2 Burnshield Dressing 600 mm x 400 mm

QTY / CASE 4



5 items

+ Belt Attachment

SUITABLE FOR
Outdoor
Travel
Mining

PORTA PACK

CODE	900805
SIZE IN MM	250 x 35 x 130
DETAIL	Nylon bag. Belt attachment. Items individually labelled with instructions for use.
CONTENTS: 5 ITEMS	x1 Burnshield Dressing 100 mm x 100 mm x1 Burnshield Dressing 200mm x 200 mm x1 FAD No. 2 (50 mm x 75 mm) x1 FAD No. 3 (75mm x 100 mm) x1 Pair of Gloves

QTY / CASE 10



14 items

PVC Case

SUITABLE FOR
Small to medium burns
Home
Travel

EASycARE KIT

CODE	550012
SIZE IN MM	200 x 63 x 130
DETAIL	PVC easy to carry plastic box. Items individually labelled with instructions for use.
CONTENTS: 14 ITEMS	x2 Burnshield Dressing 100 mm x 100 mm x1 Burnshield Dressing 200 mm x 200 mm x3 Burnshield Hydrogel Sachet 3,5ml x2 Conform Bandage 50 mm x3 Plasters Fabric x1 Plastipore Plastic Tape X1 Pair of Gloves x1 First Aid Scissors (Metal)

QTY / CASE 20



30 items

+ D Handle

SUITABLE FOR
Work
Travel
Home

ESSENTIAL BURN KIT MEDIUM

CODE	900819
SIZE IN MM	220 x 60 x 160
DETAIL	Medium sized and convenient nylon bag with D-handle, for minor burns and scalds. Items individually labelled with instructions for use.
CONTENTS: 30 ITEMS	x2 Burnshield Dressing 100mm x 100 mm x1 Burnshield Hydrogel Tube 25 ml x1 Burnshield Digit Dressing 25 mm x 500 mm x4 Burnshield Hydrogel Sachet 3,5 ml x2 Conform Bandage 50 mm x15 Waterproof Plasters x1 Plastipore Plastic Tape x1 Pair of Gloves x1 First Aid Scissors (Metal) x1 Rescue Blanket x1 Burnshield Pocket Guide

QTY / CASE 10



18 items

+ Hook Straps

SUITABLE FOR
Work
Travel
Home

ESSENTIAL BURN KIT SMALL

CODE	900818
SIZE IN MM	185 x 50 x 105
DETAIL	Small and convenient nylon bag with hook straps, for minor burns and scalds. Items individually labelled with instructions for use.
CONTENTS: 18 ITEMS	x1 Burnshield Hydrogel Tube 10ml x1 Burnshield Digit Dressing 25 mm x 500 mm x3 Burnshield Hydrogel Sachet 3,5ml x2 Conform Bandage 50 mm x8 Waterproof Plasters x1 Plastipore Plastic Tape x1 Pair of Gloves 1 Burnshield Pocket Guide

QTY / CASE 20



42 items

+ Pullout Pouches

SUITABLE FOR
Work
Travel
Home

ESSENTIAL BURN KIT LARGE

CODE	900820
SIZE IN MM	280mm x 60mm x 210mm
DETAIL	Large sized and convenient nylon bag with D-handle and pullout pouches, for burns and scalds. Items individually labelled with instructions for use.
CONTENTS: 42 ITEMS	x3 Burnshield Dressing 100 mm x 100 mm x1 Burnshield Dressing 200 mm x 200 mm x1 Burnshield Digit Dressing 25 mm x 0,5 m x1 Burnshield Limb Dressing 50 mm x 1 m x6 Burnshield Hydrogel Sachet 3,5ml x1 Burnshield Hydrogel Bottle 50ml x2 Conform Bandage 50 mm x20 Waterproof Plasters x1 Plastipore Plastic Tape x1 Cohesive Soft 1 Foam Bandage x2 Pair of Gloves x1 First Aid Scissors (Metal) x1 Rescue Blanket Adult x1 Burnshield Pocket Guide

QTY / CASE 5



SUN BURN HAPPENS

*EVEN TO THE
VIGILANT SPF
DEVOTEES*



EMERGENCY BURN CARE

Cools and soothes sun burn
Minimises trauma and skin damage
Safe for use on children
Non-toxic



BURNSHIELD®

www.burnshield.com | [@BurnshieldZA](https://www.instagram.com/BurnshieldZA) [f](https://www.facebook.com/BurnshieldZA) [y](https://www.youtube.com/BurnshieldZA) [in](https://www.linkedin.com/BurnshieldZA)