

Model name:	S-CAPEPLUS Evacuation Mattress – Basic Model
Article number: REF	5N88888886
GMDN:	66353
Basic UDI-DI:	8720726730S-CAPEPLUSZT
GTIN:	8720726730044
Classification fire retardancy:	M1 / EN597
Anti bacterial treatment:	Sanitized® TH 27-24
Size mattress (L x W x H):	210 x 60 x 7 cm / 82.6" x 23.5" x 2.8" 2,5 kg / 5.5 lbs
Size packed in carrying bag:	60 x 35 x 60 cm / 23.5" x 13.5" x 23.5" inches
Size packed in Wall Cover:	100 x 60 x 17 cm / 39.3" x 23.5" x 6.69" inches
Weight:	6,5 kg 14.5 lbs
Safe working load:	200 kg 441 lbs
Material:	100% Polyester
Laundry instructions:	Clean by hand: Use a wipe and disinfectant In the washing machine: Remove the cover and place in the washing machine. You can wash the cover with the following washing instructions: Washing temperature, 40 degrees (Celsius), dry at a low temperature, do not iron, do not dry clean, do not bleach, do not use washing powder with fabric softener
Warranty:	S-CAPEPLUS has one-year guarantee, from the invoice date, on confection and manufacturing faults.
Declaration:	Produced in alliance with the European Regulations: CE MD The product is classified in Class 1 according to Regulation MDR 2017/745/EU. Applied standards: NEN-EN-ISO 13485:2016 EN NEN-EN-ISO 14971:2019 EN

